



VantixTM Research System

User Manual

Universal Sensors Ltd

Abbey Barns
Duxford Road
Ickleton
Cambridge
CB10 1SX
UK

Tel: +44 (0) 1799 533160

Fax: +44 (0) 1799 531385

Email: info@universalsensors.co.uk

Web: www.universalsensors.co.uk

INTRODUCTION.....	2
KEY FEATURES	2
ADVANTAGES AND BENEFITS	2
APPLICATIONS	3
SENSOR FORMAT AND STRUCTURE	3
ASSAYS.....	4
POTENTIOMETRIC DETECTION.....	5
CONFIGURING THE VANTIX™ RESEARCH SYSTEM	7
IMMOBILIZATION OF BIOMOLECULES ONTO THE BIOSENSOR POLYMER LAYER	7
SPOT COATING PROCEDURE (MANUAL)	7
TROUGH COATING PROCEDURE	9
SENSOR STRIP STORAGE	10
SENSOR HANDLING	11
STRIP CUTTING.....	11
SENSOR HOLDING DURING ASSAYS.....	12
TYPICAL ASSAY PROCEDURE.....	13
ASSAY PROCEDURE STEP-BY-STEP	13
OPTIMISATION PROCESS.....	15
SANDWICH ASSAYS	15
COMPETITIVE ASSAYS	15
REVERSE COMPETITIVE ASSAYS.....	16
VARIABLES WHICH CAN BE OPTIMISED.....	17
THE VANTIX™ RESEARCH READER.....	18
VANTIX™ RESEARCH READING DEVICE	18
HOW TO TAKE A MEASUREMENT	19
VANTIX™ RESEARCH SOFTWARE.....	20
DATA OUTPUT	23
TROUBLESHOOTING.....	25
APPENDIX 1 BUFFERS AND SOLUTIONS.....	26
APPENDIX 2 VANTIX™ RESEARCH SOFTWARE INSTALLATION.....	28
APPENDIX 3 EXAMPLE ASSAY.....	33

Introduction

Biosensors are devices which detect the amount of analyte present in a solution. They consist of an immobilised biological component in contact with a transduction device that converts a biological signal into a quantifiable electrical signal. The Universal Sensors System allows for very sensitive, rapid and highly reproducible bioassays to be performed. The system uses a potentiometric biosensor with a conductive polymer as a solid phase for the attachment of the bioassay and as a sensing element. An immunocomplex is formed during the incubation step which contains an electrochemically active label. USL use a patent protected method of measurement: sensors are immersed in a buffer containing an analyte, a label generates a shift in potential either directly or inversely proportional to the amount of analyte in the sample and is measured by the reader device. Typical assays based on this system have a 5-10 minute incubation step and a two minute reading step. Assays can be customised to produce sensitivity over a wide range of analyte concentrations, from pg to µg/ml as desired.

Key Features

The Vantix™ Research system is made up of three main components:

- Sensor strips upon which the bioassays are performed.
- Reader device into which the sensors are inserted.
- Software which performs the measurement and generates the results output.

Advantages and Benefits

- Robust and reliable.
- Very short incubation times (typically 5 to 10 minutes total).
- Easy to use and simple format.
- Rapid optimisation processes.
- Many assay formats are possible, as long as they can incorporate an enzyme such as HRP to generate the electrical potential.
- Small sample volumes and amounts of reagents.
- Measurement of potential change is more sensitive than colorimetric, chemiluminescent or fluorescence measurement.
- Precise quantitative measurements in complex sample matrices e.g. blood, serum and urine, with high sensitivity and at a low cost.
- Little or no sample preparation.
- Universality and flexibility-
 - Many electrochemically active reactions can be easily incorporated either onto the sensor surface or into the required reagents.
 - Proteins, nucleic acids, toxins, drugs, contaminants, cells, bacteria and viruses are all viable analytes.
- Accurate and reproducible.
- Long term stability.
- Low cost and suitable for automation.

The Universal Sensor technology combines the advantages of both ELISA and lateral-flow immunoassays (LFIAs) without their disadvantages. It has better sensitivity and quantification than ELISA and is as robust, straightforward and as fast as LFIA. It avoids the key disadvantage of ELISA; that of extensive sample preparation and lengthy incubation times. An average ELISA will take 3 hours

to complete. A full Vantix™ - based assay can be run in less than half an hour. The Vantix™ assay also avoids the main disadvantage of LFIA; low sensitivity and a lack of quantification. Compared with other assay formats, the key advantage of Vantix™ is its ease of use. Very few preparation or handling steps are necessary and consequently it can be easily performed by unskilled users and used on-site during sample collection. The combination of simplicity and the rapidity mean that testing is extremely cost-effective.

Applications

The Vantix™ assay is suitable for the rapid, quantitative, low levels (ppt) detection of a wide range of analytes (see list below). The system is also suitable for research labs wanting to transfer existing tests or develop new tests. It can also be used for kinetic experiments of immobilized enzymes and inhibition studies.

Because of the simplicity of the system and short incubation times the assay development and/or optimisation can take just few hours when using high quality bioreagents.

Feasible analytes include:

- Proteins
- Disease markers
- Nucleic acids
- Glucose, urea, creatinine
- Drugs/metabolites
- Digoxin, antibiotics, steroids
- Cells
- Bacteria
- Virus's
- Chemical contaminants e.g. PCB, Dioxins, PAH, etc
- Heavy metals
- Ions K⁺, Mg⁺⁺, etc

Sensor Format and Structure

Sensors are manufactured using established screen-printing technology. The sensor elements are printed layer by layer starting with conducting tracks followed by reference and working electrodes. Finally a dielectric is printed on top of other layers leaving defined openings on the reference and working electrodes and insulating conducting tracks (Figure 1). For the majority of tests a conducting polymer is deposited electrochemically onto the working electrode forming the sensing element and immobilisation matrix of the working electrode.




Figure 1. USL sensor physical features. Part A of this image shows a sensor strip consisting of 40 sensors. Part B shows a close up of 3 pairs of working and reference electrodes.

Assays

Immunoassays normally make use of one or more enzymes linked to one or more of the reactants to allow quantification through change in a particular physical property: for example, the development of a colour after the addition of a suitable substrate/ chromogen in a conventional optical ELISA (Figure 2).




Figure 2. Schematic of a sandwich ELISA.

In sandwich immunoassays, the analyte is 'sandwiched' between two different antibodies. The amount of product (e.g. coloured substance) is proportional to the concentration of antigen in the sample.

In competitive assays, an analyte and a labelled analyte (competitor) compete for immobilized antibody binding sites (Figure 3). The amount of bound labelled competitor is inversely proportional to the amount of analyte in the sample. The signal (e.g. change in colour) is inversely proportional to the analyte concentration.




Figure 3. Competitive and Sandwich assays. Schematic of the complexes formed and the reaction occurring at the surface of the biosensor for both competitive and sandwich assay formats.

Potentiometric Detection

In Vantix™ assays, the electrical potential difference between the 2 electrodes is measured, rather than a colour change. The difference in potential is created by the reaction which takes the place on or near the surface of the working electrode. One of the many processes taking place near the working electrode is the transfer of electrons directly from the electrochemically reactive electrode to the substrate molecule via the active site of the enzyme. (Figure 4). The signal is a result of several processes, which take place during the enzymatic reaction.




Figure 4. The reaction catalysed by Horse Radish Peroxidase (HRP). OPD converting into DAP (coloured product in optical tests).

A sensor with a formed immunocomplex is placed into the measuring substrate solution containing an H-donor, for example OPD (o-phenylenediamine), with hydrogen peroxide as the active substrate. The products generated by the enzymatic reaction change the electrochemical properties of the polymer layer, which in turn shifts the resting potential of the polymer layer (patent protected method of measurement). The difference between the two potentials is related to the concentration of analyte in a sample. The more HRP bound to the complex the greater the shift from resting potential.

Configuring the Vantix™ Research System

Immobilization of Biomolecules onto the biosensor polymer layer

The first step in the assay complex formation on the sensor is immobilising the first biomolecule to the polymer coated working electrode via physical adsorption. Rapid immobilisation protocols that work well for medium to large proteins have been developed. Smaller proteins may need to be conjugated to other molecules e.g. BSA or biotin (and incubated with sensors coated with streptavidin). Post absorption treatment with commercial stabilising buffers, drying and correct storage allow sensors to be stored for months to years if required.

Sensors are provided in strips of 40, which are ready to coat (Figure 1). Close examination of the working electrode (black lower electrode) shows that there is a 1.0 mm circular gap in the dielectric layer directly over the electrode. This is the area that is coated with the bioreceptor(s). The other electrode (reference) should not come into contact with the first immobilised molecule(s).

There are two approaches to coating the sensor strips, either the immersion of the working electrode end of the strips into liquid containing the first biomolecule held in a trough or microtitre plate well, or the spotting of the biomolecules on to the circular gap of the working electrode.

Spot Coating Procedure (manual)

1. Prepare a 20-50 µg/ml stock solution of the first biomolecule in coating buffer. We recommend 50mM potassium phosphate buffer pH 8.0 (appendix 1) for IgG and molecules with isoelectric points of 7 or above and 50mM citrate buffer pH 5.0 for molecules with isoelectric points below 7 (appendix 1).
2. Lay the sensor strips on a rigid piece of cardboard (electrode gap sides facing upwards) and fix them with a small piece of tape attached to the ends of the strip.
3. Using a low volume pipette carefully spot 3 µl drops of the biomolecule stock directly onto the 1 mm working electrode gap (note the drop must cover the gap entirely, Figures 5 and 6).
4. Without disturbing the droplets, move the strips on the cardboard into a 37 °C incubator for 15 minutes (note the droplets should not be permitted to dry out during this time, 3 µl is normally sufficient not to evaporate over 15 minutes).
5. As soon as the 15 minutes have elapsed, wash the working electrode end of the strips in a trough containing approximately 50ml of 50mM Potassium Phosphate buffer pH 8.0 by dipping them in and wiggling/swirling several times. It is OK to expose both electrodes to the washing solution, but the contact pads should remain dry. If they get wet they should be wiped dry. Wipe the sensors on tissue to remove buffer. Repeat 3 times.
6. Immerse the working electrode end of the strips into a trough of stabilising buffer for 10 minutes (use USL stabiliser or commercially available Stabliguard Choice Stabiliser, refer to Appendix 1 for ordering details).
Note: use a volume of stabilising buffer that allows only the working electrode to be fully covered by the solution, but not the reference electrode (Figure 7)
7. Remove from trough, wipe the back of the strips (but not the front with the working electrodes on) and fix onto the cardboard (electrode gap side facing upwards).
8. Dry for 2 hours in an incubator at 37 °C before prolonged storage, or for 30 minutes at 37 °C before immediate use. Note that the spotting part of this process can be automated using a suitable liquid handling robot e.g. a BioDot machine. Alternatively USL can provide sensors spotted with biomolecule and stabilised to your specification.




Figure 5. Spot coating. A sensor strip being manually coated with a 3 μ l drop of biomolecule.




Figure 6. Spot coating close up. The 3 μ l drop is placed directly over the 1.0 mm gap of the working electrode.



Figure 7. Sensors in stabilising buffer. Note that only the black working electrode is immersed in the stabilising buffer, not the reference electrode (buffer has been coloured yellow for ease of viewing).

Trough Coating Procedure

Biomolecules can be immobilised by immersing the working electrodes into the coating solution containing biomolecule(s). We recommend using uncut strips, but it is possible to use pre-cut sensors inserted into a holder as shown on Figure 8. This method results in greater uniformity of coating but has higher reagent consumption, which is partly compensated by using lower concentrations and longer incubation times.

Suggested immobilisation concentrations are 1-20 µg/ml, typically equivalent to microtitre plate based ELISA.

1. Fill a trough with enough biomolecule coating solution to cover the black working electrode but **NOT** the reference (silver) electrode (Figure 8). The buffers used in this method are the same as those suggested in the manual spot coating section.
2. Incubate the strips in this solution for 15-30 minutes at room temperature.
3. Wash the working electrodes end of the strips in a trough containing approximately 50 ml of 50mM Potassium Phosphate buffer pH 8.0 by dipping them in and wiggling/swirling several times. It is OK to expose both electrodes to the washing solution, but the contact pads should remain dry. If they get wet they should be wiped dry. Wipe the sensors on tissue to remove buffer. Repeat 3 times.
4. Immerse the working electrode end of the strips into a trough of stabilising buffer for 10 minutes (use USL stabiliser or commercially available Stabliguard Choice Stabiliser, refer to Appendix 1 for ordering details).

Note: use a volume of stabilising buffer that allows only the working electrode to be fully covered by the solution, but not the reference electrode (Figure 7)

5. Remove from trough, wipe the back of the strips (but not the front with the working electrode) and fix onto the cardboard (electrode gap side facing upwards).
6. Dry for 2 hours in an incubator at 37 °C before prolonged storage, or for 30 minutes at 37 °C before immediate use. Note that the spotting part of this process can be automated using a suitable liquid handling robot e.g. a BioDot machine. Alternatively USL can provide sensors spotted with biomolecule and stabilised to your specification.




Figure 8. Trough coating. The sensor working electrodes are immersed in a trough containing biomolecule.

Sensor Strip Storage


Once sensors are coated, dried for 2 hours and sealed (e.g. in an opaque bag with desiccant), they can be stored at 4 °C for extended periods of time. The time of possible storage will depend on storage conditions and biomolecules stability. For example antibodies have been shown to retain their function for several years when stored in optimal conditions (basically without access of light, moisture and oxygen).

We recommend that the strips are kept stuck onto paper bases in a re-sealable plastic bag inside a sealed foil bag containing desiccant and that these are stored at 4 °C until use. This will protect from light, temperature and moisture. Bags can be opened to remove sufficient sensors for an assay, then resealed and returned to 4 °C for future use of the remaining sensors.


Sensor Handling

Strip Cutting.


Prior to use in an assay (and post coating) sensors should be cut from the strips (Figure 9). Any remaining unused sensors can be returned to the packaging for future use.



i) Trimming the base.



ii) Trimming the silver electrodes.



iii) Cutting individual sensors.

When coated sensors are required for an assay, the following method is the best way for preparing sensors for use. Use a good quality, sharp pair of scissors.

- Trim the base so that about 3 mm of plastic remains below the working electrode (i).
- Cut off the end excess plastic so that the end sensors have about a 2-3 mm gap from the silver track to the end.
- Carefully trim the top of the strip so that the silver tracks are cut evenly at the level of the top of the silver track of the reference electrode – cut the whole strip (ii).
- Finally cut out individual sensors by cutting between a reference electrode silver tracks of one sensor and the working electrode silver tracks of the neighbouring sensor (iii).

The above should yield individual sensors as shown (iii). Any sensors that have been cut to within 1 mm or less to the electrodes or silver tracks should be discarded i.e. these components must be surrounded by a minimum of 1 mm of plastic.

Note it is important to ensure the top part of the sensor is trimmed as shown.

Figure 9. Preparation of individual sensors from strips. The process of trimming sensors is described above.

Sensor holding during assays

Simple sensor holders are provided, which help to simplify sensor manipulation during the assay steps. The contact ends of the sensors are inserted into the narrow holes in the holder. Note it is good practice to insert the sensors in the correct orientation (as shown), with the working electrode on the left side and the reference electrode on the right side of the sensors (Figure 10).




Figure 10. Sensors holder. The left image shows the contact end of a sensor being inserted into one of the narrow holes of the holder. The right hand image shows the holder with 12 sensors, note the correct orientation of the sensors. The spacing of the slits in the holder matches the spacing of wells in standard SBS 96-well microtitre plates. Fewer than 12 sensors can be used if desired.

The holder allows for the easy manipulation of the sensors. They can easily be inserted into assay plates and washed between incubation steps (Figure 11).



Figure 11. Sensors in a holder inserted into 12 wells of an assay plate.

The working electrodes are immersed into the wells of an assay plate, containing the required bioreagents, for incubation steps. We recommend the use of flat bottom 96 well polypropylene assay plates as this is a low protein binding plastic, which is critical for assays requiring sensitivity below 1 ng/ml. In a typical assay, sensors are immersed in either 200 or 300 μ l volumes of reactive mixtures for incubation steps and washed in troughs between the steps.

Note: 200 µl covers the working electrode while 300 µl covers both the working and reference electrodes. It is important that the reference electrode is either completely immersed or completely dry.

Incubation times vary between 5-15 minutes depending on the system, but typically 5-10 minutes of incubation is sufficient. Note it is important to use exactly the same incubation times for all inter-assay work (unless optimising assay timings), as the signals will vary with varying incubation lengths resulting in potentially misleading results.

Typical assay procedure


It has been mentioned above that a typical assay development and basic optimisation can be achieved within just a few days if good quality reagents are used. The process of assay development based on Vantix™ sensors is quite different from the one applied for conventional ELISA. Because of short incubation times the whole assay run can be done within 40 minutes including a bioreceptor coating step. Therefore there is no need to set up many variable conditions at the starting point as is the standard practice when developing ELISAs. We recommend starting with just 3-6 sensors and one or at the very most two reagent dilutions and taking the results into consideration in further development steps. Often people are tempted to “check” other conditions straight away, which is fine. However, it actually saves time, reagents and effort to develop the assays by performing consecutive “whole” runs and narrowing down conditions. The final decision on development is of course ultimately down to the particular researcher. We recommend starting with defining the target time for the assay, e.g. 15 minutes and reconsider the time when/if all other development options have been exhausted.

Assay procedure step-by-step

1. Make up reading solution and set up the reader apparatus.
2. Cut out previously coated sensors and insert each into a sensor holder.
3. Prepare the samples in assay buffer; we would typically start with few samples with a wide dilution range (e.g. a zero sample and the top of your expected analyte concentration range).
4. Place the samples (200 or 300 µl volumes) into wells of the same row of a flat bottom 96 well plate.
5. Insert the sensors into the samples (ensure that the working electrode is entirely immersed) and incubate for 5 -10 minutes. While sensors are incubating, prepare 2 troughs with approximately 40 ml assay buffer and tissues laid out (to wash the sensors with at the end of incubation).
6. Start the reading software and enter any parameters.
7. When the incubation time has elapsed, immediately wash the sensors by removing the sensors from the assay plate, and briefly swirling the electrode ends in the trough of assay buffer and patting (blotting) both sides of the sensors dry on tissue (Figure 12).
8. Repeat the wash/blot step three more times.
9. Repeat the washing steps in the second trough of assay buffer.
Note: During washing both electrodes and the part above them up to connecting pads can be immersed and washed. The washing step is critical for most assays.
10. If an additional incubation step is required (e.g. sandwich assays) place the sensors strip in the assay plate row with the reagent and incubate for 5-10 minutes, then repeat the full washing process using fresh assay buffer.
11. After the end of the final incubation & wash step, immediately proceed to the reading step.




i) Fill two troughs with 40-50 ml of wash buffer and layout several low lint tissues such as Kimberly-Clark medical wipes.




ii) Remove the holder with the sensors from the assay plate, insert the electrodes into the wash buffer to approximately halfway up the sensors and swirl for several seconds.

Note it is important to perform this wash then blot step first rather than blot then wash.



iii) Gently wipe the back side of the sensors on the tissue paper.



iv) Gently blot the front side of the sensors on the tissue, repeat this process twice using the first trough of buffer, then for an additional 3 times using the second trough of buffer. It is best to blot/wipe the sensors in a fresh part of the tissue each time.

Figure 12. Washing procedures for sensors fixed in a holder.

Optimisation process

Like all assay systems, assays using the Vantix™ Research system will require optimisation. The advantages this system is that the assays are very quick to run, so assays can be optimised in a relatively short time. The objective is to maximise the signal range between the extreme values of analyte concentration you wish to measure. Such difference can be between zero and the lowest concentration (LOD) or cut-off value (the point at which the analyte is exceeding the maximum allowed level). Some clinical assays require a stable calibration curve within a certain range of concentrations, e.g. clinical window.

It is very important to identify the targets and to know exactly what needs to be achieved BEFORE starting the assay development.

The signals mentioned below are using the USL reading solution (see appendix 1), other reading solutions will give different signals.

Optimisation should always be carried out using the matrix you intend to use in the final assay, or solution of similar properties. An assay developed in a different matrix could behave very differently.

Note: there are many formats for sandwich, competitive and sequential assays available for consideration. If you have little experience in assay development it is recommended that you refer to immunoassay literature to study the subject prior to starting experiments.

For the Vantix™ Research system incubation times should not be reduced below 5 minutes. Manual steps and deviations in incubation times will contribute to standard deviations (SDs).

Also keep in mind that the Vantix™ Research system relies on diffusion exchange between the electrode and bioreagent, which takes at least two minutes to complete.

Sandwich assays

Capture antibodies should initially be immobilised at a concentration of 20-50 µg/ml (the exact concentration can be optimised once the subsequent steps have been optimised). Initially use three analyte concentrations – zero, the desired minimum to be detected and the expected maximum. Incubate with standards (known analyte preparations) for 15 minutes and wash. Incubate with the second antibody (conjugated to HRP) for 5-10 minutes at a starting concentration of 1 µg/ml (or 1/1000 dilution if the concentration is not known), wash and read. If this result is not satisfactory 1) vary the conjugate concentration, 2) change the incubation times, 3) review coating concentrations, 4) review buffers/solutions. Once this initial assay is working, additional standard concentrations should be incorporated to continue optimisation.

Ideally in the OPD based USL measuring solution “zero” should be at -120 mV +/-20 mV. The separation between zero and the lowest concentration (LOD) should be at least 15 to 20 mV, the signals for the highest analyte concentration should be from +40 to +50 mV.

Competitive assays

Capture antibodies are initially immobilised at a concentration of 20-50 µg/ml (this concentration can be optimised once the subsequent steps have been optimised). We do not recommend using coating concentrations below 10 µg/ml as this can affect the calibration curve over prolonged storage due to deterioration of antibodies on the surface. Using good stabilising solution can help reduce this problem significantly.

Start with defining the incubation times. In competitive assays it is not always beneficial to have longer incubations. We recommend starting with 10 minutes for a one step competitive assays. We recommend 7 minutes for the first incubation and 5 minutes for the second incubation in sequential assays.

Start with zero, low and high concentrations of the analyte in assay diluent buffer or a specific matrix. Keep the concentration of the competing analyte (usually the analyte conjugated to HRP) constant. The starting concentration of the competitor should be approximately in the centre of the desired analyte concentration range. Mix reactants, add the sensors and incubate as directed, then wash the sensors and read.

Ideally in USL OPD based reading solution the signal for zero analyte should be from 20 to 50 mV and for the highest concentration of analyte from -80 to -100 mV.

Optimise the assay by 1) changing conjugate concentration, 2) changing incubation times, 3) changing coating concentration, 4) considering alternative formats.

Reverse competitive assays

Some reagents work best in this format where the antigen (or antigen-conjugate) is immobilised on the sensor with the antibody and the sample being incubated with the sensor. As in the standard competitive assay format described above, the concentration of analyte is inversely proportional to the signal gained. The antibody can be conjugated to HRP or a second incubation with an HRP conjugated antibody raised against the animal species of the first antibody can be used for detection. An initial coating concentration of 20-50 µg/ml for the analyte conjugated to an inert protein (such as BSA) should be used, and approximately 0.1 µg/ml for the antibody in solution should be used. The final concentrations will need to be optimised for each reagent set.

Competitive assay			Sandwich assay		
Analyte (unit/ml)	mV value	Δ mV	Analyte (unit/ml)	mV value	Δ mV
0	21.0	0.0	0	-87.5	0.0
10	-8.2	29.3	10	-65.2	22.
50	-47.0	68.1	50	-47.0	40.4
100	-65.2	86.3	100	-8.2	79.3
300	-87.5	108.6	300	21.0	108.6

Table 1. Examples of acceptable signal values (mV) for the two assay formats. Units can be pg/ml, ng/ml, µg/ml, ppt, ppb or ppm, depending on assay set up and bioreagent specifications. For competitive assays the signal gained decrease as the analyte concentration increases. With sandwich assays the signal gained increases as analyte concentration increases. The standard USL measuring solution gives a measuring range from -120 mV to 50 mV, while other proprietary solutions give other results, e.g. Arista TMB has a range from -60 mV to over 100 mV. Note that although the range given by commercial substrate solutions may appear superior to that of USL OPD based measuring solution, commercial solutions tend to give much higher intra and inter assay variations. USL is currently developing a proprietary measuring solution which will improve on the abilities of the OPD-based measuring solution. This will be released soon but samples are available on request. Some ready-to-use substrate solutions are not compatible with the Vantix™ system, mainly because they contain solvents and electrochemically active additives. Please consult USL prior to the use of any non-recommended substrate.

Variables which can be optimised

- Each of the biomolecules involved e.g. capture antibody concentration.
- Conjugate concentration.
- Incubation times.
- Incubation temperature.
- Presences of blocks in buffer e.g. BSA, Casein.

Details of optimisation strategies are covered in the Universal Sensors document 'Instructive Training Toolkit'.

The Vantix™ Research Reader

Vantix™ Research Reading Device

The Vantix™ Research system allows up to 12 sensors to be read at the same time. The contact pad end of each sensor is inserted into the base of the reader head (Figure 15), both electrodes of the sensors are then immersed in the measuring solution contained in a trough (Figure 14), note it is important that both electrodes are fully covered by the solution. The electrical potential is measured via the Vantix™ software on a PC or PDA device for 1 to 2 minutes.




Figure 13. Vantix™ Research reader system. The figure shows a measurement being taken on the Vantix™ software, collecting data from 12 sensors in the reading device. The system allows the electrical potential of up to 12 sensors to be read at the same time. Measurements take typically 1 to 2 minutes. Real time readings and graphs are generated during the read and results are exported in CSV files.




Figure 14. Vantix™ Research measuring device. Twelve sensors inserted into the Vantix™ device.




Figure 15. The head of Vantix™ Research device. 12 sensors inserted in the 12 reading slots. Note the correct orientation of the sensors and that they are inserted to the same depth.

How to take a measurement

(Please refer to “Vantix™ Research software” section for more detailed information)

Once the assay has been completed and the sensors washed for the final time, they should be individually removed from the holder, inserted into the measuring device, immersed into the measuring solution and read. It is a good practice to keep all steps reproducible in terms of times and actions. The first two steps below should be completed prior to the end of the assay so that the sensors can be measured immediately after the final wash.

1. Open the software and choose from the reading options (the system will automatically open with the previously saved protocol). Modify the assay setup as required and press “run the program” to go to the next window.
2. Make up and pour 50 ml of reading solution into a 50 ml trough (or 10 ml into a 10 ml trough) place the trough as shown above (Figure 14).
3. Wipe the sensors dry after the last washing step (it will not affect the sensors). Remove the first sensor from the holder (ensure that the contact ends are dry, if not dry carefully with tissue paper), insert the contact end into the first position of the Vantix™ reading head. Each sensor must be positioned correctly and inserted firmly into the reading head (push the sensor carefully into the holder until it can go no further). Repeat with the remaining sensors, keeping them in the same order and maintaining the correct orientation (refer to Figure 15).
4. Insert the legs of the reading device head into the reader base holes so that both of the sensor electrodes on each of the sensors are fully immersed in the measuring solution (Figure 14). Ensure that the sensors are not touching the base or sides of the trough (adjust the legs where required).
5. Press “start run” in the live data window/ control board box (see Figure 17). At the end of the reading, a dialogue box appears prompting the user to name and save the read file.
6. The measuring solution can be reused up to five times over a 4 hour period. However, it is a good practice to use a fresh portion of solution by the end of assay optimisation and/or while accumulating statistical data. Protect it from light and evaporation between readings e.g. surround with foil.

Vantix™ Research Software.

The software allows the user to record the potential generated by sensors in reading solution. Recording happens in real time and can be recorded over a set time period with set frequency. Data can be collected has raw signal data or as relative data compared to a zero control in the first channel. Signals can be collect in any number of channels from 1 to 12 and holding potentials can be applied during recording if required. The software can also be used to apply a calibration curve to data and output sensor measurements as concentrations. Finally, once parameters have been set they can be saved as an .apr file for instant recall.

Note: The Vantix™ Research device should be plugged in and the software started **before** placing sensors into the Vantix™ Research device. A small current is generated in the Vantix™ Research device when the software is started and this can potentially damage any sensors plugged into the Vantix™ Research device.

The first window generated when opening the Vantix™ Research software is shown in Figure 16. It allows the user to set up the reading parameters.




Figure 16. Vantix™ Research software first window. Shows the first window generated when opening the Vantix™ Research software.

Settings: Allows the user to choose the method by which data is displayed. “**Absolute values**” will record the raw data of signals measured at control point 1, control point 2 and on the last second (length) of the experiment. Make sure all three time points (control point 1, control point 2 and length) differ by at least one second (The standards for these are control point 1 = 60s, control point 2 = 90s and length = 91s or 120s). The “**Relative values**” function will display values zeroed against the first channel for control point 1 and control point 2. The “Relative value” function is often used in competitive assays with a “zero sensor” (incubated in a negative control sample) inserted into channel 1. The values from other channels are then reported relative to a negative control.

Sensor selection and identification: The user can select the reader channels used (deselecting channels that will not contain sensors) and name the channels if desired. If no name is given to a channel the default naming is Sensor 01 through 12.

Potential step: This feature allows the user to set a holding voltage on the sensors for a set time (normally 10-15 seconds). In some applications it helps to minimise background signals prior to taking a reading. If this feature is used it is recommended that the holding voltage is set at a value equal to the mV signal given by an antibody coated sensor (e.g. in a sandwich assay with **no** conjugate/reporting antibody). The typical value in USL reading solution is -115 to -150 mV. To disable this feature the user should insert a '0' (zero) into the “Num. of steps” box.

Note: if using a potential step we highly recommend also using the “current limit” feature, set to between 20 and 100 nA. More detailed instructions for this will be given below.

Measurement parameters: The user defines the time of the reading by choosing “length” of measurement, typically 60-120 seconds. The time between each reading sample can also be defined. In the example above readings will be taken every second (1000 milliseconds). The result for each second is reported in a CSV file generated by the software. Two additional results can be displayed in the “live” window and subsequently in the reporting table of the exported CSV file. The control point 1 and control point 2 allow the user to define at what time the software takes an additional reading. The time for each control point, and the length of the read must be different by at least one second, e.g. length = 121s, control point 1 = 60s, control point 2 = 120s.

In most cases the user will not need require the control point functions.

Calibration tables: This allows individual dual calibration curves for each sensor, allowing multi-analyte testing. This feature is only essential for an assays being prepared for commercialisation. Currently it is being validated by USL, and has not been fully optimised in this version of software.

Stored parameters: Once parameters have been set this feature can be used to save and recall them for repeat use. At USL we call it a “protocol” for a specific assay or repeated runs with the same parameters.

Run the program: When this button is pressed two other windows will be displayed. In order to return to the first window use the “review parameters” button. After any changes are made the user must press the “run the program” button again to set the changes.




Figure 17. VantixTM Research software second window. Shows the windows that open once the 'run the program' button is clicked in the first reader window.

Graph window: During a reading this displays the real-time millivolt signals generated for each sensor against time.

Graph re-scaling sliders: These allow the user to modify the range and scale of the graph during a run.

Control point window: This displays the millivolt readings for the control points.

Live data window and control board The ADC mV shows the real-time millivolt reading value for each sensor, the adjacent colour matches the line colour displayed in the graph window during the run (and can be changed by clicking the word colour and selecting a new one). By pressing the 'start run' button the program begins recording the run as it has been set in the first window. The 'start run' button should only be clicked once the sensors are loaded into the reading head and the reading head placed into the

base (with the sensor electrodes immersed in reading solution). The 'stop run' button can be used to prematurely stop the run (the data will not be saved).

The current limit option should be ticked and set between 20 and 100 mA if a potential step is being used. Sensors may be damaged and/or give higher variations if a potential step is applied without setting a current limit.

Data output

The Vantix™ Research software generates a CSV file at the end of the reading step. The user is prompted to name and save the file in a location of their choice. An example of the output is shown below.

Absolute values

File name:	H:\Duncan.Smith\ASSAY\Expt25.csv				
System time:	12/02/2008 15:36				
Name of the experiment:	Expt25				
Sensor name	Raw Data, 121s	Control point 60s	Result 1	Comment	Control point 120s
Sensor 01	-0.56	-0.23			-0.55
Sensor 02	-16.19	-15.83			-16.17
Sensor 03	-36.94	-38.2			-36.92
Sensor 04	-75.69	-80.2			-75.7
Sensor 05	-109.75	-110.42			-109.74
Sensor 06	-106.44	-107.39			-106.42
Sensor 07	-49.97	-52.39			-49.99
Sensor 08	-49.84	-51.89			-49.86
Sensor 09	-50.41	-52.36			-50.39
Sensor 10	-53.22	-54.77			-53.2
Sensor 11	-57.59	-61.11			-57.61
Sensor 12	-54.72	-57.86			-54.74

Relative values

Sensor name	Raw Data, 121s	Control point 60s	Result 1	Comment	Control point 120s
Sensor 01	-0.56				
Sensor 02	-16.19	15.63			15.62
Sensor 03	-36.94	36.38			36.37
Sensor 04	-75.69	75.13			75.15
Sensor 05	-109.75	109.19			109.19
Sensor 06	-106.44	105.88			105.87
Sensor 07	-49.97	49.41			49.44
Sensor 08	-49.84	49.85			49.84
Sensor 09	-50.41	49.85			49.84
Sensor 10	-53.22	52.66			52.65
Sensor 11	-57.59	57.03			57.06
Sensor 12	-54.72	54.16			54.19

Table 2 and 3. Vantix™ Research software data output. The tables below show parts of an output file (exported as a CSV file and opened using Microsoft excel). The top table shows output using the absolute values function. The lower table shows output using the relative values function. In each case the top two rows show the location of the file and the time it was generated. The third row shows the name of the file. Column 1 contains the names of each of the sensors. Column 2 contains the final millivolt measurement for each sensor at the end of the measuring time. Columns 3 and 6 show the control point readings taken during the run e.g. at 60 and 120 seconds (note: these times are defined by the user in the first software window (Figure 16)).

Not shown in the tables is the additional data collected e.g. for every second of the reading for each sensor (note: the recording period is defined by the user in the first software window (Figure 16)). This data can be used for further detailed examination of the result.

The columns with “result” and “comment” are generated by the user in calibration tables in “settings”. It can be used for comments such as “fail”, “pass”, etc. in commercial assays.

Troubleshooting

	Observation	Explanation	Resolution
1	Very high or very low (about +/-1000mV) values in live window.	The electrodes are not covered by the reading solution, e.g. ref. electrode is out.	Re-adjust the “legs” of the reading head so that both electrodes are properly immersed in the solution.
		Sensor not inserted correctly into reading head.	Stop read, remove and reinsert the sensor and then restart read.
		Software was opened after immersing sensors into the reading solution.	The sensors must be discarded and the assay re-run with fresh sensors (refer to software handling above).
2	All or some sensors read zero.	The silver electrodes at the top of the sensor is not cut correctly.	Cut the silver electrode to produce two separate electrodes (refer the Fig. 9).
		There is a strip of contact pad from the neighbouring sensor.	Cut out the strip of the neighbouring contact pad.
3	Wavy graph during reading.	Reference electrode not fully immersed in reading solution.	Top up reading solution or/and adjust the height of the reading head leg nuts.
4	High backgrounds are observed in assay.	May be due to non-specific binding of conjugate.	<ol style="list-style-type: none"> 1. Titrate conjugate until acceptable background values are observed. 2. Reduce incubation times and or 3. Add block e.g. 1% BSA to assay buffer.
5	Weak response is observed in assay.	Zero analyte in competitive assays should result in high signals (e.g. greater than +20 mV for USL reading solution). If signals are lower than this then the concentration of conjugate may be too low.	<ol style="list-style-type: none"> 1. Increase the concentration of conjugate. 2. Increase the incubation time. 3. Immobilised analyte may not have bound efficiently, the active site may have immobilised preferentially to the surface or the conjugate may have denatured – refer to immobilisation strategy (page 8).
6	Poor passive adsorption/immobilisation resulting in low response.	Some reagents will not readily immobilise to the surface or the binding site may bind to the surface effectively blocking it.	<ol style="list-style-type: none"> 1. Biotinylate the reagent and incubate with 80ug/ml streptavidin coated sensors. 2. Conjugate the reagent to an other protein e.g. BSA. ?? 3. If the reagent is an antibody, immobilise an anti-species antibody to the sensor surface, then incubate with the reagent. 4. Use a buffer close to the isoelectric point of the molecule. 5. Increase the concentration of the molecule being immobilised.
7	Dynamic range of assay suboptimal (the gap between response for zero and response for very high concentration is small).	Assay may require further optimisation.	<ol style="list-style-type: none"> 1. Optimise further each of the key components e.g. immobilised reagent, concentration of conjugate. 2. Increase the assay time. 3. Try block solutions (e.g. 1% BSA) - may reduce non-specific binding. 4. Consider the ratio of HRP to antigen used in the manufacture of the conjugate, increasing the HRP to antigen ratio may improve sensitivity. 5. As with ELISAs there is a maximum sensitivity due to the biochemistry of the reagents; consider trying alternative reagents.
8	Poor duplicates (the values for the same concentration are far apart).	Manual spotting of the first reagent can give elevated CVs.	<ol style="list-style-type: none"> 1. Try the trough coating method and compare results. 2. Check that the droplets don't dry out during the 15 minute coating step at +37 °C. 3. Proper washing of the electrodes is very critical, review the washing protocol. 4. Keep the timing as strict as possible, e.g. add components using multichannel pipette, keep the time between the steps constant, etc.

Table 4. Troubleshooting. Problems that may be encountered with the Vantix™ Research system and suggested resolutions.

Appendix 1 Buffers and solutions

Assay buffer and Coating buffer

The coating buffer is 50 mM Potassium Phosphate, pH 8.0.

The assay buffer is 0.1 M Potassium Phosphate, pH 7.8.

From the premade stock solutions: 0.5 M Potassium Phosphate monobasic and 0.5 M Potassium Phosphate dibasic, dilute the required portions of both buffers in separate containers to obtain the required molarity (e.g. 1/10 for 50 mM). Then add monobasic buffer to dibasic buffer stepwise with mixing until the required pH is obtained.

Stabilising solution for proteins on the surface

USL surface stabilising solution.

Alternatively, use Surmodics product: StabilGuard Choice, Biomolecule stabiliser, Cat# SG-02.

Stabilising solution for proteins/antibodies and conjugates in solution

Normally, 1-3% BSA or casein in Potassium Phosphate buffer or PBS is good enough for R&D experiments.

Alternatively, use Surmodics product: Stabilzyme-HRP Cat# SZ-02.

USL Reading solution

(Reading values are typically from -130 mV to +40 mV, but can be from -150 mV to +60 mV).

Add 2.5 ml of Sodium Perborate (formulation given below) solution to 50 ml of basic solution (formulation given below) and mix well. Store refrigerated in a tin foil covered container (solution is light sensitive).

Preparation of the USL reading solution constituent components:

50 mM Sodium Citrate buffer pH 5.0

Prepare stock solutions of 0.2 M Sodium Citrate and 0.2 M Citric Acid and mix them to gain a solution of pH 4.85. Mix one part of this with three parts MilliQ water to get 50 mM Citrate buffer with pH 5.0. Store at 4 °C for up to 1 month.

Sodium Perborate solution

Dissolve 120 mg of Sodium Perborate in 20 ml of 50 mM Citrate Buffer pH 5.0.

Store refrigerated in a tin foil covered container (light sensitive). Make fresh on the day of measurement.

Basic solution

Dissolve one 20 mg OPD tablet in 200 ml 50 mM Citrate buffer, pH 5.0

Store in a tin foil covered container (light sensitive) at RT. Make fresh on the day of measurement.

Other reading solutions

Sigma Fast Reading solution

(Reading values are typically from -60 mV to +80 mV, but can be from -80 mV to +120 mV).

Sigma fast tablets can be used. Make up to 50 ml in milliQ water rather than the recommended 20 ml (Cat# P-9187).

TMB-based reading solution

(Reading values are typically from -50 mV to +100 mV, but can be from -60 mV to 130 mV).
We recommend Arista Biological TMB (Cat# ARTMB-0010).

Note: Generally third party reading solutions are optimised for use in ELISAs, they are less stable than USL reading solution and must be changed more frequently, e.g. after two measurements. You can use a smaller 10ml trough to minimise the volume used.

Blocks

(Normally not required).

Some reagents are sticky to the sensors and give non-specific binding. The addition of 2% BSA or 1-2% casein to the assay buffer is usually sufficient to eliminate non-specific binding. In some cases a 5 minute pre-incubation of the sensors (bearing the first reagent) in block before incubation with the analyte can be beneficial. Alternatively the block can be added to the coating or stabilising solutions as well as to the assays diluents/solutions.

Reagents/components

OPD tablets – Sigma, P-5412.

Na-Perborate – Sigma (Aldrich), 372862.

Citric Acid – Sigma, C-7129.

Na-Citrate – Sigma, S-4641.

K-Phosphate dibasic – Sigma, P-3786.

K-Phosphate monobasic – Sigma, P-0662.

BSA- Sigma, A7906.

Casein, www.sdt-reagents.de, CBC1.

SigmaFast tablets – Sigma, P-9187.

TMB – Arista, ARTMB-0010.

Biomolecule stabiliser for proteins on the solid phase, StabilGuard Choice, Surmodics (Diatec AG for Europe), Cat#SG-02.

Stabilising solution for HRP-conjugates, Stabilzyme-HRP, Surmodics (Diatec AG for Europe), Cat#SZ-02.

96 well plates – Appleton Woods, CC713.

50ml troughs – Appleton Woods, BIO201.

Appendix 2 Vantix™ Research Software Installation

Please follow carefully!

If you have used the Vantix™ system on your computer previously and think you may be running a previous version of the software just replace the “051212 uts MFC” file from V-UTS folder on the C: drive and desktop of your computer with the “nUTS 080924” file on the installation CD.

If you have never used Vantix™ system on your computer follow the installation guide below.

Software Installation

The Vantix™ Research system from Universal Sensors Ltd is designed to measure the responses of sensors produced by Universal Sensors Ltd. The Vantix™ unit connects to a PC through a standard USB connector.

Before use the software must be installed on the host computer from the supplied CD.

Computer Requirements

The programs for operating the Vantix™ Research system have been developed for PC's running Microsoft Windows XP. The programs may work on Windows 2000 or later, however use on operating systems prior to Windows 2000 cannot be guaranteed.

Installation of Drivers

The Vantix™ Research system uses a FT232BM chip to communicate with the host computer. The drivers for this chip need to be installed before the Vantix™ Research system can communicate with the host computer.

System Components

Prior to connecting the Vantix™ Research system to a PC, ensure you have the following components:

- Vantix™ Research sensor head
- Vantix™ Research base plate
- USB cable
- Drivers & Software CD


Do not connect the Vantix™ Research sensor head to the computer until instructed to do so!

Copying required files to the PC


Note: The following screen captures are only guides, and different screens may be seen on your PC during the installation due to your system configuration.

- Create a folder 'V-UTS' in the root directory of the drive (normally C:) on the required PC e.g. C:\V-UTS.
- Copy all the files from the supplied CD into the new folder.
- Connect the smaller plug of the USB cable into Vantix™ Research sensor head.
- Now connect the larger USB plug into a spare USB port of the PC. This will launch the Windows 'Found New Hardware' Wizard. This may appear either as a speech bubble at the bottom of the screen, or the system may immediately display a dialog box as shown below. If only the speech bubble is visible, click on the bubble to show the dialog box.

- If the dialog box is as shown in the picture below, click on the 'No, not this time' option and then click 'Next'.




- Select the 'Install from a list or specific location (Advanced)' option as shown below and then click "Next".




- Select 'Search for the best driver in these locations'. If currently checked, deselect the 'Search removable media (floppy, CD-ROM...)' check box. Click on the 'Include this location in the search:' and enter the file path in the combo-box ("C:\V-UTS\Drivers\R9052154") or browse to


this directory by clicking the browse button and then selecting the R9052154 folder. Once the file path has been entered, click 'Next' to proceed.




- If Windows XP is configured to warn when unsigned drivers are about to be installed, a screen may be displayed informing the user that the drivers have not passed the Windows Logo testing scheme. Click on 'Continue Anyway' to continue.
- A screen will be displayed as Windows XP copies the required driver files. After the files are copied, the following screen will be seen. Click 'Finish' to end the wizard.




- Windows XP will now install the COM port emulation files. The 'Found New Hardware' Wizard will begin again. Click on the 'No, not this time' option again and click 'Next' to continue.



- Select the 'Install from a list or specific location (Advanced)' option as shown below and then click "Next".



- As previous, select 'Search for the best driver in these locations'. The previously entered options and file path should be shown. If so click 'Next' to continue otherwise enter the file path in the combo-box ("C:\V-UTS\Drivers\R9052154") then click 'Next'.



- Windows XP will copy the files required for the serial port emulator. The following dialog box will eventually be displayed. Click 'Finish' to end the wizard.



- Make a copy of "nUTS 080924" and paste on Desktop (do not use a shortcut!).
- Start the program by double-clicking on the shortcut or the icon of the program

Appendix 3 Example assay

Mouse IgG Assay

This is an example of the results gained from a typical sandwich assay using well characterised bioreagents. Sensors were coated with goat anti mouse IgG (Fc specific) Ab at 20 ug/ml in 50 mM potassium phosphate buffer (pH 8.0). Mouse serum was diluted in 100 mM potassium phosphate buffer (pH 7.8) +1% BSA to give IgG concentrations between 0.1 and 1000000 ng/ml. Sensors were incubated in 300 µl samples of each serum dilution for 5 minutes then washed in 100 mM potassium phosphate buffer (pH 7.8). Sensors were then incubated in goat anti mouse IgG HRP conjugate diluted 200x in 100 mM potassium phosphate buffer (pH 7.8) +1% BSA and washed again. Sensors were then plugged into the Vantix™ Research head and immersed in universal sensors reading solution. Signals were recorded for 2 minutes and then plotted against a log₁₀ IgG concentration scale.




Figure 18. A mouse IgG standard curve. This standard curve was gained using typical sandwich assay conditions on the Vantix™ Research system.

The results shown in Figure 18 show a linear range between 10 and 10000 ng/ml. This result was gained following just two 5 minute incubations and a 2 minute read. The details of the optimisation can be found in another USL document entitled 'Instructive training toolkit'.